

## INFORMATION ASSAYS HEMATO-ONCOLOGY LAB - UNIT MOLECULAR DIAGNOSTICS EMC

Abnormality	Percentage Malignant Cells*	VAF#	Assay
<i>RUNX1-RUNX1T1</i>	1%		RQ-PCR1
<i>CBFB-MYH11</i>	10%		RQ-PCR1/RT-PCR
<i>FLT3</i> ITD	10%	0.05 (ratio)^	PCR > LMA
<i>FLT3</i> TKD	10%	0.05 (ratio)^	PCR > LMA
<i>NPM1</i> mutation	0.01% (RQ-PCR <sup>‡</sup> )	1% (NGS)	RQ-PCR1/NGS1
<i>CEBPA</i> mutation	10%	5%	NGS2
<i>ASXL1</i> mutation	10%	5%	NGS1
<i>TP53</i> mutation	10%	5%	NGS1/NGS3
<i>RUNX1</i> mutation	10%	5%	NGS1
<i>KIT</i> mutation	10%	5%	NGS1
<i>IDH1</i> mutation	10%	5%	NGS1
<i>IDH2</i> mutation	10%	5%	NGS1
<i>PML-RARA</i>	1%		RT-PCR
<i>MLL-AF4</i>	1%		RT-PCR
<i>BCR-ABL</i> (breakpoint)	1%		RT-PCR
<i>BCR-ABL</i> (e13/e14-a2)	0.001%		RQ-PCR2
<i>BCR-ABL</i> mutation	only >0.01%		Sanger sequencing
<i>NUP214-ABL</i>	0.1%		RQ-PCR1
<i>FIP1L1-PDGFR</i>	1%		RT-PCR
<i>SETBP1</i> mutation	10%	5%	NGS1
<i>NRAS</i> mutation	10%	5%	NGS1
<i>CSF3R</i> mutation	10%	5%	NGS1
<i>BRAF</i> mutation	10%	5%	NGS1
<i>SF3B1</i> mutation	10%	5%	NGS1
<i>MYD88</i> mutation	10%	5%	NGS1
<i>JAK2</i> V617F mutation	1%	2%	RQ-PCR1/NGS4
<i>JAK2</i> exon12 mutation	10%	5%	NGS4
<i>MPL</i> mutation	10%	5%	NGS4
<i>CALR</i> mutation	10%	5%	NGS4
<i>NRAS</i> mutation	10%	5%	NGS1
<i>EZH2</i> mutation	10%	5%	NGS1
<i>U2AF1</i> mutation	10%	5%	NGS1
<i>CBL</i> mutation	10%	5%	NGS1
<i>SRSF2</i> mutation	10%	5%	NGS1
<i>TET2</i> mutation	10%	5%	NGS1
<i>UBA1</i> mutation	10%	5%	NGS5

\*De result not detectable in case of AML is only reliable above the indicated percentage of malignant cells or VAF  
 #VAF: variant allele frequency (eg. VAF 50% means that all the cells carry a heterozygous mutation) – VAF below the percentage indicated are also reported  
 ^The *FLT3* ITD and *FLT3* TKD mutation results are reported as ratio (*FLT3* ITD/ *FLT3* wild type of *FLT3* TKD/ *FLT3* wild type) according to ELN 2017 [no VAF!]  
 ‡NPM1 mutant type A,B and D

PCR: polymerase chain reaction  
 RT-PCR: reverse transcriptase – PCR  
 RQ-PCR1: real-time quantitative – PCR (in-house test)  
 RQ-PCR2: real-time quantitative – PCR (Cepheid)  
 LMA: length mutation analyse  
 NGS1: Next generation sequencing Illumina Trusight Myeloid panel  
 NGS2: Next generation sequencing custom panel *CEBPA*  
 NGS3: Next generation sequencing custom panel *TP53*  
 NGS4: Next generation sequencing custom panel *JAK2, CALR, MPL*  
 NGS4: Next generation sequencing custom panel *UBA1*

On the next page the markers determined in the Hemato-oncology laboratory unit Molecular Diagnostics are indicated per malignancy.

Disclaimer next generation sequencing:  
 Neutral polymorphisms, silent mutations, likely not-pathogenic variants and variants, that are currently unclassifiable, will not be reported. The method is limited for the detection and/or exclusion of mutations below a variant allele frequency of 5% as well as (larger) deletions, insertions and indels. Some parts of genes are more difficult to sequence such that mutations may be missed.

Markers Hemato-oncology laboratory unit Molecular Diagnostics ranked by malignancy

<b>AML and MDS-EB</b>	<b>Type mutations/ region</b>
<i>RUNX1-RUNX1T1</i>	4bp insertie/ exon12
<i>CBFB-MYH11</i>	complete gene
<i>NPM1</i> mutation	exon14 en exon15
<i>CEBPA</i> mutation	D835/ exon20
<i>FLT3-ITD</i>	exon12
<i>FLT3-TKD</i>	complete gene
<i>ASXL1</i> mutation	exon2-11
<i>RUNX1</i> mutation	exon2,8-11,13,17
<i>TP53</i> mutation	R132/ exon4
<i>KIT</i> mutation	R140 en R172/ exon4
<i>IDH1</i> mutation	
<i>IDH2</i> mutation	
<b>CLL</b>	
<i>TP53</i> mutation	exon2-11
<b>CMML</b>	
<i>ASXL1</i> mutation	exon12
<i>SETBP1</i> mutation	(part of) exon4
<i>NRAS</i> mutation	codons 12,13 en 61/ exon2 en exon3
<i>RUNX1</i> mutation	gehele gen
<b>CNL and atypical CML</b>	
<i>ASXL1</i> mutation	exon12
<i>CSF3R</i> mutation	exon14-17
<i>SETBP1</i> mutation	(part of) exon4
<b>Hairy cell leukemia</b>	
<i>BRAF</i> mutation	V600/ exon15
<b>Mastocytosis</b>	
<i>KIT</i> mutation	exon2,8-11,13,17
<b>MPN</b>	
<i>CALR</i> mutation	exon9
<i>JAK2</i> mutation	V617F en exon12
<i>MPL</i> mutation	W515/ exon10
<b>M. Waldenström</b>	
<i>MYD88</i> mutation	L265/ exon3-5
<b>MDS</b>	
<i>SF3B1</i> mutation	K700/ exon13-16
<i>EZH2</i> mutation	complete gene
<i>U2AF1</i> mutation	exon2 en exon6
<i>ASXL1</i> mutation	exon12
<i>RUNX1</i> mutation	complete gene
<i>TP53</i> mutation	exon2-11
<i>CBL</i> mutation	exon8 en exon9
<b>MDS-RS</b>	
<i>SF3B1</i> mutation	K700/ exon13-16
<b>VEXAS</b>	
<i>UBA1</i> mutation	exon3

Predispositie genpanel (uitsluitend na overleg)

informatie: <https://hema13.erasmusmc.nl/diagnostiek.html>

Disclaimer next generation sequencing:

Neutral polymorphisms, silent mutations, likely not-pathogenic variants and variants, that are currently unclassifiable, will not be reported. The method is limited for the detection and/or exclusion of mutations below a variant allele frequency of 5% as well as (larger) deletions, insertions and indels. Some parts of genes are more difficult to sequence such that mutations may be missed.

### Assay description

The variants in the genes indicated are determined by next generation sequencing using the Illumina TruSight Myeloid panel (TSCA design (probe hybridisatie, extensie-ligatie en PCR)) or custom panels (nested PCR) on the Illumina MiSeq<sup>1</sup>. The coverage of all target genes is at least 500x. Of the exon-intron boundaries 5bp is analysed. The variants are identified with an in-house developed software pipeline. In this software read alignment is done with BMAP<sup>2</sup> and SAMtools<sup>3</sup>, whereas 'variant calling' with MuTect<sup>4</sup>, GATK<sup>5</sup>, Varscan<sup>6</sup>, Indelocator<sup>7</sup> and Pindel<sup>5</sup>. All analyses are carried out on two independent DNA samples. Only the DNA variants present in both samples are reported. To classify the variants as pathogenic several *in silico* programs, such as SIFT<sup>8</sup>, Polyphen<sup>9</sup> en Mutation Taster<sup>10</sup>, Franklin [https://franklin.genoox.com] en de Genome Aggregation Database (gnomAD) are used.

#### Genomic regions with poorer coverage

chr1:115256418-115256622  
chr1:36932117-36932335  
chr21:36259341-36259561  
chr3:128200640-128200861  
chr3:128204781-128204979

#### Part of the following gene

NRAS  
CSF3R  
RUNX1  
GATA2  
GATA2

### References

1. <https://emea.illumina.com/>.
2. BMAP short-read aligner, and other bioinformatics tools (<http://sourceforge.net/projects/bbmap/>).
3. Li H, Handsaker B, Wysoker A, et al. The Sequence Alignment/Map format and SAMtools. *Bioinformatics* 2009;25:2078-9.
4. Cibulskis K, Lawrence MS, Carter SL, et al. Sensitive detection of somatic point mutations in impure and heterogeneous cancer samples. *Nat Biotechnol* 2013;31:213-9.
5. DePristo MA, Banks E, Poplin R, et al. A framework for variation discovery and genotyping using next-generation DNA sequencing data. *Nat Genet* 2011;43:491-8.
6. Koboldt DC, Zhang Q, Larson DE, et al. VarScan 2: somatic mutation and copy number alteration discovery in cancer by exome sequencing. *Genome Res* 2012;22:568-76.
7. Indelocator. at <http://www.broadinstitute.org/cancer/cga/indelocator>.
8. Ng PC, Henikoff S. SIFT: Predicting amino acid changes that affect protein function. *Nucleic Acids Res.* 2003;31(13):3812-4.
9. Adzhubei I, Jordan DM, Sunyaev SR. Predicting functional effect of human missense mutations using PolyPhen-2. *Curr Protoc Hum Genet.* 2013;Chapter 7:Unit7.20.
10. <http://www.mutationtaster.org/>.

#### Disclaimer next generation sequencing:

Neutral polymorphisms, silent mutations, likely not-pathogenic variants and variants, that are currently unclassifiable, will not be reported. The method is limited for the detection and/or exclusion of mutations below a variant allele frequency of 5% as well as (larger) deletions, insertions and indels. Some parts of genes are more difficult to sequence such that mutations may be missed.